Infection control quote worksheet

Facility Name:	 Date:
Facility Contact:	
Arjo Representative:	

Please check all equipment that you would like included in your quote.

Flites									
Flite Style — Clip									
	Product Name	Product Descriptions	Product Code	Size	Quantity				
O	General Purpose Clip Flite for Arjo passive lifts	Disposable clip flites, designed to be used for one patient during their stay. Specify Size : S, M, L, XL, XXL. 10 per box.	MFA1000M-X-L1						
O	Sara 3000® and Sara Flex Clip Flite	Single patient use <i>Sara 3000</i> clip flites combines infection control and safe patient handling. 10 per box	MFA3000-L1	One Size					
O	Clip Limb Flite	Single patient use clip limb flite assists the caregiver with lifting, supporting, and positioning limbs. (47" x 5.5") - 5 per box	VIG22000-L1	One Size					
O	Turning/Repositioning Flite	Single patient use turning/repositioning flite assists the caregiver for supine turns in bed. (31.5" x 17.7" x 5.5") - 5 per box	VIG22003-L1	One Size					
Flite Style -	– Loop								
O	General Purpose Loop Flite for Arjo passive lifts	Single patient use disposable loop flite supports infection control and optimizes comfort. Specify Size : S, M, L, XL, XXL. 10 per box.	MFA2100-X-L1						
O	Loop Limb Flite	Single patient use loop limb sling assists the caregiver with lifting, supporting, and positioning limbs (47" x 5.5") - 5 per box	VIG22001-L1	One Size					
O	Bariatric Repositioning Sling	Single patient use bariatric repositioning sling for supine position transfers. (145 cm x 210 cm, 57" x 82") - 5 per box	VIG22004-L1	One Size					
•	Repositioning Sling	Single patient use disposable sling for repositioning and moving the patient up the bed. (208 cm x 119 cm, 82" x 47") - 5 per box	AHD001-L1	One Size					
O	Sara Plus® Flite	Single patient use <i>Sara Plus</i> Flites combine infection control and safe patient handling routines. 10 per box	MFA4000-L1	One Size					
		a 3000/ Limb Flite Turning / Repositioning Fli	General Purpose te Loop Flite	Repositioning S	Sling				





Infection control quote worksheet

Please check all equipment that you would like included in your quote.

Wipeable Standing Clip Slings								
Our Wipeable Standing Clip Sling is nonporous and has welded seams. The body, strap, and attachment points are designed to be able to be wiped clean between use, helping to reduce infection risks associated with shared patient handling devices. Avoids the need for laundering. Suitable for use with Sara Flex and Sara 3000.								
	Product Name	Product Description	Product Code	Size	Quantity			
O	Wipeable Standing Clip Sling for Sara 3000 and Sara Flex	Standing clip sling for optimal infection control between patients. Suitable for general sit to stand, transfers, and toileting. Specify Size : S, M, L, XL.	TSS.50XW-L1					
Wipedown Standing Slings								
The Wipedown sling is a two piece system with no velcro, helping to reduce infection risks associated with shared patient handling devices. Avoids the need for laundering. Suitable for use with Sara Plus.								
0	Wipedown Sling for Sara Plus, Encore, Chorus	Wipe down sling for optimal infection control between patients. Suitable for general sit to stand, transfers, and toileting. Specify Size : S, M, L, XL.	KKA5420-X-L1					
		Wipeable Wipedown Standing Clip Sling						

If you have any questions please contact your Arjo sales representative or if completed return to customer service via fax at (888) 594 2756 or email us.cc@arjo.com

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure injuries and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.





